NEW EDITION 2017

HP Special

Swiss Dental Diamond
since 1968 by Vanetti SA
www.vanetti.ch


HP Special


Utensili diamantati HP Instruments diamantés HP Diamant - Schleifinstrumente HP
Diamond instruments HP Instrumentos diamantados HP Алмазные инструменты HP

ISO 104

PM - HP - H - PM - PM


106-125 µm
HP
STANDARD
ISO 104

REF.-No	G801	G830	G859	G368
ISO 1/10 mm	065	085	065	085
Fig.-No ISO
Head L mm	13,5	21,5	23,5	22
Prices P	14	16	15	15

	801	830	859	368
	065	085	065	085

	13,5	21,5	23,5	22

Imballaggio • Emballage • Verpackung
Package • Embalaje • Упаковка

1

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⌚ Vitesse recommandée 2000-5000 min⁻¹

(GB) User information's

After the resin has become hard, the rims have to be reduced with the usual rotary instruments at an angle of 90 degrees with regard to the inside of the baseplate or the prostheses. With the help of the usual rotary instruments you need to touch up till you get the required thickness. The new diamond instrument DIA-TESSIN offers a real solution to this problem. It is also used at an angle of 90 degrees with regard to the inside of the baseplate or the prostheses and you get rounded rims with the right dimension without any difficulties. This step-by-step method can be applied both to the upper and lower jaw. The new DIA-TESSIN rotary instrument can be used with any resin; it is also suitable for rim-correction of the prostheses, for preparation of relinings, etc.

Summarised, the basic advantages of the new instrument DIA-TESSIN are:

- A conic head
- Two grooves with specific dimensions: al 2,2 mm, bl 3,2 mm.

These two grooves allow to respect the prescribed dimensions with great accuracy. It is also very useful for the dentist who wants to make some corrections to the baseplate or the prostheses without any loss of dimension and with always rounded rims.

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